



Hospital Chemical Cleaner List

Cleaners

Chemical Name	Manufacturer
Ammonia (15% solution)	Generic
CideZyme	ASP J&J
CideZyme XTRA	ASP J&J
Clean Screen Wipes	Clorox
Easy Screen Cleaning Wipe	PDI
EmPower, EmPower Foam	Metrex
Formula 409 cleaner Degreaser Disinfectant	Clorox
Formula 409 Glass & Surface Cleaner	Clorox
Formula 409 Heavy Duty Degreaser	Clorox
GreenWorks Cleaners	Clorox
GreenWorks Manual Pot & Pan Detergent	Clorox
MetriClean2	Metrex
MetriLube	Metrex
Metrex	Metrex
MetriShine	Metrex
MetriWash	Metrex
MetriZyme, MetriSponge	Generic
Mild dishwashing detergent	Clorox
Pine-Sol	Ecolab
Sekusept Cleaner	Clorox
Tilex Instant Mildew Remover	Clorox
Tilex Soap Scum Remover & Disinfectant	Clorox
Urine Remover	Metrex
VioNexus Foaming soap w/Vitamin E	Generic

Disinfectants

These chemicals are used to reduce contamination and the spread of bacteria and germs.

Product	Manufacturer
1-Stroke Environ	Steris
4 in one Disinfectant & Sanitizer	Clorox
Accel INTERvention RTU & Wipes	Diversey
Accel PREvention Concentrate	Diversey
Accel PREvention Wipes & RTU	Diversey
Accel TB	Diversey
Alpha HP	Diversey



Avert	Diversey
Avert Sporidical Disinfectant Cleaner Wipes	Diversey
Birex Quat™ Wipes	Biotrol
BirexSE	Biotrol
Bleach (<50% solution)	Generic
Bleach (10% solution)	Clorox
Bleach (Germicidal-Hypochlorite)	Generic
Bleach Germicidal Disinfectants, Wipe & Liquid	Clorox
Broad Spectrum Quaternary Disinfectant Cleaner	Clorox
Cavicide 1, CaviWipes 1	Metrex
Cavicide Bleach, CaviWipes Bleach	Metrex
Cavicide, CaviWipes	Metrex
CDC-10	Spartan Chemical
Cidex	ASP J&J
Cidex OPA	ASP J&J
Citrace Hospital Disinfectant and Deoderizer	Clorox
Clean Quick Broad Ranger Quaternary Sanitizer	P&G
CleanCide RTU	Wexford Labs
CleanCide Wipes	Wexford Labs
Clean-Up Disinfectant Cleaner w/Bleach	Clorox
Clinell Universal Spray	Clinell
Clinell Universal wipes	Clinell
Comet Disinfecting - Sanitizing Bathroom Cleaner	P&G
Comet Disinfecting Cleaner w/Bleach	P&G
Consume Eco-Lyzer	Spartan Chemical
Control III® Disinfectant and Germicide	Maril Inc.
Control III® Laboratory Germicide	Maril Inc.
Coverage Plus Germicidal Wipes	Steris
Coverage Plus NPD Disinfectant	Steris
Coverage Spray HB+ Disinfectant	Steris
Coverage Spray TB Disinfectant Cleaner	Steris
Coverage Spray TB+ Disinfectant Cleaner	Steris
Diffense	Spartan Chemical
Disinfecting BioStain & Odor Remover	Clorox
Disinfecting Wipe	Clorox
Dispatch Hospital Cleaner Disinfectant Towells w/Bleach	Clorox
End Bac II	Diversey
Envirocide	Metrex
Environ LpH se Disinfectant	Steris
Environ LpH st Sterile Disinfectant	Steris
Environ Vesphene se Disinfectant	Steris
Environ Vesphene st Disinfectant	Steris
Envy Foaming Disinfectant Cleaner	Diversey
Envy Liquid Disinfectant Cleaner	Diversey



Expose II 256	Diversey
EZ-Kill Quat Alcohol Cleaner Disinfectant Wipes	Clorox
Fuzion Clean & Disinfectant	Clorox
Germicidal Bleach	Clorox
Glutaraldehyde (3.4% solution)	Generic
GS HD Disinfectant 256	Spartan Chemical
GS Neutral Disinfectant Cleaner	Spartan Chemical
H2O2 cleaner disinfectant	Generic
Halt	Spartan Chemical
Hepacide Quat II	Spartan Chemical
Hydrogen Peroxide Cleaner Disinfectants, Wipe & Liquid	Clorox
Hydrogen Peroxide Disinfecting Cleaners	Clorox
Incides N	Ecolab
Incidin	Ecolab
Incidin Active	Ecolab
Incidin Foam	Ecolab
Incidin Liquid	Ecolab
Incidin Oxywipe/Oxyfoam	Ecolab
Incidin Oxywipe/Oxyfoam S	Ecolab
Incidin Plus	Ecolab
Incidin Pro	Ecolab
Incidin Rapid	Ecolab
LpH ag Disinfectant	Steris
LpH III se Phenolic Disinfectant	Steris
LpH III st Phenolic Disinfectant	Steris
LpH se One Step Disinfectant	Steris
Metricide	Metrex
Metricide OPA Plus	Metrex
Metricide Plus 28	Metrex
Metricide Plus 30	Metrex
Metriguard	Metrex
Micro-Cide	Micro Scientific
Mr. Clean Professional Disinfecting Multi-Purpose Cleaner	P&G
Multi-Surface Quat Alcohol Cleaner Disinfectant Wipes	Clorox
Neutral Disinfectant Cleaner	Ecolab
Neutral Quat Disinfectant Cleaner	3M
Opti-Cide Max	Micro Scientific
Opti-Cide3	Micro Scientific
Opti-Cide3 RI	Micro Scientific
Oxivir 1 RTU	Diversey
Oxivir 1 wipes	Diversey
Oxivir Five 16	Diversey
Oxivir Plus	Diversey
Oxivir Tb RTU	Diversey



Oxivir Tb wipes	Diversey
Oxycide Daily Disinfectant and Cleaner	Ecolab
P&G Pro Line Disinfecting Floor & Surface Cleaner II	P&G
P-4D	Spartan Chemical
PerCept	Diversey
PerDiem	Diversey
PreSept	ASP J&J
Process LpH st Sterile Disinfectant	Steris
Process NPD One Step Disinfectant	Steris
Process NPD st Sterile Disinfectant	Steris
Process QDS Disinfectant	Steris
Process Vesphene list Disinfectant	Steris
Quat alcohol cleaner	Generic
Quaternary Disinfectant Cleaner	Ecolab
Rescue RTU Gel & Wipes	Diversy
RTU Wex-Cide	Wexford Labs
Sani-Cloth Active	Ecolab
Sani-Cloth AF3 Germicidal Disposable Wipe	PDI
Sani-Cloth Bleach Germicidal Disposable Wipe	PDI
Sani-Cloth HB Germicidal Disposable Wipe	PDI
Sani-Cloth Plus Germicidal Disposable Cloth	PDI
Sani-Cloth Prime Germicidal Disposable Wipe	PDI
Sani-Prime Germicidal Spray	PDI
Sani-T-10	Spartan Chemical
Sekusept AKTIV	Ecolab
Sekusept Extra N	Ecolab
Sekusept Multienzyme P	Ecolab
Sekusept Pulver Classic	Ecolab
Sekusept Pure Clean	Ecolab
SekuseptT Plus	Ecolab
Spartan Hard Surface Disinfecting Wipes	Spartan Chemical
Spic and Span Disinfecting All-Purpose Spray & Glass Cleaner	P&G
Spor-Klenz RTU Cold Sterilant	Steris
Sterigent	Spartan Chemical
StorOx 2.0	P&G
Super Neutral L	Spartan Chemical
Super-Sani cloth Germicidal Disposable Wipe	PDI
XXX TB-Cide Quat	Spartan Chemical
Thymo-Cide	Wexford Labs
Triad III	Diversey
VersaSure Cleaner Disinfectant Wipe	Clorox
Vesphene II SE* One Step Disinfectant	Steris
Vesphene III se Phenolic Disinfectant	Steris
Vesphene III st Phenolic Disinfectant	Steris



Vesta Syde SQ Disinfectant	Steris
Vesta Syde SQ st Disinfectant	Steris
Virex II 256	Diversey
Virex Plus	Diversey
Virex TB RTU	Diversey
Virox AHP 5 RTU	Diversey
Wex-Cide 128	Wexford Labs
Wex-Cide RTU	Wexford Labs
Wexford Disinfectants Wipes	Wexford Labs
Wonder Woman Wipes (IPA)	Emtec

Hand Cleaners

Hand cleaners, while often times low in abrasiveness and having only mild astringent qualities, can still be problematic with the plastics used in medical devices and should be considered in testing prototypes and choosing materials.

Product	Manufacturer
Healthy Soap	Purell
Opti-Scrub	Micro Scientific
Wex-Lan Skin Cleanser	Wexford Labs

Hand Lotions

Product	Manufacturer
Hand Lotion	Eucerin
HealthLink AloeSooth Moisterizing Lotion	Clorox
Sunscreen - Sport Performance SPF50	Banna Boat
VioNex Hand lotion	Metrex

Hand Sanitizers

Product	Manufacturer
Advanced Hand Sanitizer	Purell
Alcare Foamed Handrub	Steris
Alcare OR Foamed Handrub	Steris
Alcare Plus Foamed Handrub	Steris
Alcohol Hand Wipes	Steris
Cal Stat Plus Handrub	Steris
HealthLink AloeGuard Antimicrobial Soap	Clorox
HealthLink GBC AloeGel Instant Hand Sanitizer	Clorox
Kindest Kare Foam Handrub	Steris
Sani-Hands	PDI



Skinman Clear	Ecolab
Skinman Complete	Ecolab
Skinman Complete Pure	Ecolab
Skinman Foam	Ecolab
Skinman Soft	Ecolab
Soft 'N Sure Gel	Steris
Soft 'N Sure Handrub	Steris
SpitAcid	Ecolab
VioNex Antimicrobial Liquid Soap	Metrex
VioNexus Antimicrobial Foaming Soap	Metrex
VioNexus No Rinse Spray	Metrex

Household Solutions

Product	Manufacturer
Bug Spray Deep Woods Off VIII	Johnson
Castile Soap	Dr Bronner's
D Towel lotion (Hypochlorite)	unknown
Hydrogen peroxide (3% solution)	Generic
Isopropyl alcohol (70%-91% solution)	Generic
Ketchup	Heinz
Mustard	French's
Nail Polish Remover (acetone)	Generic
Nail Polish Remover (non-acetone)	Generic
Safflower Oil	Generic